



LIBRARIES
UNIVERSITY OF WISCONSIN - MADISON

Dexon advertisement.

[s.l.]: [s.n.], 1975

<https://digital.library.wisc.edu/1711.dl/UZQZTSEOCWY578B>

<http://rightsstatements.org/vocab/InC/1.0/>

The libraries provide public access to a wide range of material, including online exhibits, digitized collections, archival finding aids, our catalog, online articles, and a growing range of materials in many media.

When possible, we provide rights information in catalog records, finding aids, and other metadata that accompanies collections or items. However, it is always the user's obligation to evaluate copyright and rights issues in light of their own use.

DEXON® polyglycolic acid suture

Complete Product Information

DESCRIPTION DEXON® polyglycolic acid Suture—Synthetic, Absorbable is a homopolymer of glycolic acid. The sutures are sterile, inert, noncollagenous, nonantigenic, nonpyrogenic, flexible, braided, either undyed with a natural beige color or dyed green with D&C Green No. 6. They are uniform in size and tensile strength, but are smaller in diameter than other Absorbable Surgical Sutures of equivalent tensile strength.

ACTIONS When DEXON is placed in tissues a minimal tissue reaction occurs, which is followed by a microscopic layer of fibrous connective tissue which grows into the suture material. Animal studies (subcutaneous tissue in rats, rabbits and dogs) revealed minimal absorption at 7 to 15 days, significant absorption at 30 days, and maximum resorption after 60 to 90 days. Tensile strength, not being a function of the absorption rate, may vary from tissue to tissue, depending in part on the rate of hydrolysis. The early tensile strength of DEXON is reported to be greater than that of comparable chromic catgut. In animal studies (subcutaneous tissue of rats, rabbits and dogs) there occurred a 30 to 50% loss of initial tensile strength in one week and 60 to 80% loss in two weeks, at which time as in the case with chromic catgut no useful residual tensile strength remains.

INDICATIONS DEXON sutures are indicated whenever absorbable sutures and ligatures are employed, provided that extended approximation of tissues is not required or where undue postoperative tension or stress is not anticipated. DEXON is also indicated where non-absorbable sutures are used for their characteristic lack of tissue reaction.

CONTRAINDICATIONS DEXON is contraindicated where extended approximation of tissues under strain must be maintained.

WARNINGS The safe use of this suture in cardiovascular and neural surgery has not been established.

Under certain circumstances, notably orthopedic procedures, immobilization by external support may be employed at the discretion of the surgeon.

Do not resterilize. Discard opened, unused sutures.

PRECAUTIONS. Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds. DEXON suture knots must be properly placed to be secure. Therefore, place first throw in precise position for final knot, using a double loop, tie second throw square using horizontal tension; additional throws may be used as desired. Skin sutures which remain in place for periods of longer than seven days may cause localized topical irritation and the extended portion of suture may be snipped off after five to seven days, as indicated.

ADVERSE REACTIONS Those reactions that have been reported include tissue reaction or inflammation, fibrous or granulation tissue, wound separation and bleeding and accumulation of fluid around subcuticular stitches.

DOSAGE AND ADMINISTRATION Use as required. Do not resterilize. Discard opened, unused sutures.

HOW SUPPLIED Suture sizes 7-0 through 2. Supplied in cut lengths or ligating reels, non-needled or affixed to the various DAVIS + GECK ATRAUMATIC® needles in one and three dozen packages.

Rev. 1/75

DG
DAVIS+GECK

American Cyanamid Company
Pearl River, N.Y. 10965