



LIBRARIES
UNIVERSITY OF WISCONSIN - MADISON

Mark7 advertisement.


[s.l.]: [s.n.], 1975

<https://digital.library.wisc.edu/1711.dl/EGYCJCKV4HUFS8G>

<http://rightsstatements.org/vocab/InC/1.0/>

The libraries provide public access to a wide range of material, including online exhibits, digitized collections, archival finding aids, our catalog, online articles, and a growing range of materials in many media.

When possible, we provide rights information in catalog records, finding aids, and other metadata that accompanies collections or items. However, it is always the user's obligation to evaluate copyright and rights issues in light of their own use.



flexible, accurate, easy-to-read, disposable intrauterine sound

- **small diameter**

Same diameter as the Cu-7® intrauterine copper contraceptive — only 3 mm. to facilitate sounding of even nulliparas with little patient distress.

- **As flexible as the inserter with Cu-7 inside.**

- **accurate and legible**

After insertion of sound to the fundus, push the slide to the cervix, lift the proximal end of the slide to “lock” it in place on the sound; remove the sound and read the precise depth in centimeters at the pointer in the window of the slide.

See your Searle representative for a demonstration. It will be worth your while to place your first order with your surgical supply house through a representative, who will make you a one-time introductory offer. Mark-7 is available in single-use units, packaged sterile in 12s and 30s.

New **Mark-7**TM

brand of disposable intrauterine sound

SEARLE

Searle Laboratories

Division of G.D. Searle & Co.

Box 5110, Chicago, Illinois 60680