

Polyaction Knee Orthosis advertisement.

[s.l.]: [s.n.], 1984

<https://digital.library.wisc.edu/1711.dl/SWBHCQEJJPIWG9D>

<http://rightsstatements.org/vocab/InC/1.0/>

The libraries provide public access to a wide range of material, including online exhibits, digitized collections, archival finding aids, our catalog, online articles, and a growing range of materials in many media.

When possible, we provide rights information in catalog records, finding aids, and other metadata that accompanies collections or items. However, it is always the user's obligation to evaluate copyright and rights issues in light of their own use.

POLYACTION™ KNEE ORTHOSIS

Medial Lateral Anterior Control

Lightweight Thermoplastic

Unique Thigh Design

- Reduces slippage
- Increases medio-lateral stability
- Adjustable circumference

Durable Stainless Steel Polycentric Knee Joint

Adjustable Molded Tibial Section

Adjustability allows total contact, increased comfort and support

Patent
Pending



Two Control Options

1) Scott-Stop*

Adjustable Flexion
Extension Control



* Patent
Pending

Extension from 140° to 180°
Flexion from 70° to 190°

2) Non-adjustable

10° and 20° extension stops

Modified Polyaction™ Knee Orthosis has also been successfully used in controlling genu varum, valgus, and recurvatum

**For Further Production Information Contact Your Local Orthotist or
Call Toll Free 1-800-821-5795. (In Colorado 423-8510)**

Stock Measuring and Size Information
(Custom Fabrication Available From Cast)

Size	Calf Circumference	Thigh Circumference
Small	13 - 13 3/4"	17 - 18"
Medium	14 - 14 3/4"	19 - 20"
Large	15 - 15 3/4"	21 - 22"
Extra Large	16 - 16 3/4"	23 - 24"

NOTE: Measure Calf Circumference at largest point and thigh circumference 7" proximal to knee center.

SCOTT
ORTHOTIC LABS

POST OFFICE BOX 12207
5300 WEST 44TH AVENUE
DENVER, COLORADO 80212