



## Hollister Umbilical Cord Ligation advertisement.

[s.l.]: [s.n.], 1972

<https://digital.library.wisc.edu/1711.dl/I5OICHVNW37NO9E>

<http://rightsstatements.org/vocab/InC/1.0/>

The libraries provide public access to a wide range of material, including online exhibits, digitized collections, archival finding aids, our catalog, online articles, and a growing range of materials in many media.

When possible, we provide rights information in catalog records, finding aids, and other metadata that accompanies collections or items. However, it is always the user's obligation to evaluate copyright and rights issues in light of their own use.

at  
your  
fingertips...

secure  
umbilical cord  
ligation



When it's time to ligate the umbilical cord, a Hollister Double-Grip™ Cord-Clamp should be within reach. Its contoured finger-grips and wide jaw angle make one-hand application easy.

Hollister's Cord-Clamp has other benefits too: a hinge guard to keep even a large cord within the sealing area; firm-holding Double-Grip jaws to prevent slipping; a constant, even pressure to eliminate the dangers of seepage; and no need for belly bands or dressings. The clamp has a permanent, blind closure. When it's ready for removal—usually after 24 hours—the clamp is simply cut through at the hinge. Hollister provides the clipper.

This disposable, lightweight Hollister Cord-Clamp may be auto-claved, or it can be purchased in individual pre-sterilized packets. Write for samples and literature, on hospital or professional letter-head, please.

HOLLISTER

