



LIBRARIES
UNIVERSITY OF WISCONSIN - MADISON

Zimmer Orthopaedic Ltd advertisement.

[s.l.]: [s.n.], 1948

<https://digital.library.wisc.edu/1711.dl/KSTZHWZKAKOIA8N>

<http://rightsstatements.org/vocab/InC/1.0/>

The libraries provide public access to a wide range of material, including online exhibits, digitized collections, archival finding aids, our catalog, online articles, and a growing range of materials in many media.

When possible, we provide rights information in catalog records, finding aids, and other metadata that accompanies collections or items. However, it is always the user's obligation to evaluate copyright and rights issues in light of their own use.

Zimmer Orthopaedic Ltd

REGISTERED TRADE MARK





Manufacturers of Fracture Equipment

The TUPMAN OSTEOTOMY PLATE

Designed to facilitate Displacement Osteotomy of the Hip, the TUPMAN Osteotomy represents a significant advance in modern treatment

- ★ *It permits precise fitting and accurate adjustment of the McKee nail*
- ★ *The danger of bone fragments being fixed in distraction is overcome*
- ★ *Four angles of plate available*
- ★ *Reinforcing crimp gives strength and rigidity*



No. 525 SMO Stainless Steel (B.S.S. EN 58J) TUPMAN Plate

- | | | |
|----------|---------------|--|
| No. 525a | .. | with bracket 85° |
| No. 525b | .. | with bracket 90° |
| No. 525c | .. | with bracket 95° |
| No. 525d | .. | with bracket 100° |
| No. 428 | McKee's Nails | 3", 3 1/4", 3 1/2", 3 3/4", 4", 4 1/4" |
| No. 428a | .. | 4 1/2", 4 3/4", 5", 5 1/4", 5 1/2", 5 3/4", 6" |

Each nail is supplied with SMO locking nut, threaded 7 mm. x 1 mm. pitch.


BRIDGEND · GLAM · G.B.

Phone: BRIDGEND 938/1208 Grams: ZIMORPEDIC, BRIDGEND


London Showrooms: Zimmer House, 176-8 Brompton Road, S.W.3
KNightsbridge 1919/9672

Zimmer Orthopaedic Ltd

REGISTERED TRADE MARK



Manufacturers of Fracture Equipment




The WAINWRIGHT OSTEOTOMY PLATE

- ★ One piece solid forged stainless steel to B.S.S. EN 58J (American equivalent type 316 or SMO)
- ★ Can be used with ordinary S.P. nail impactor
- ★ Note two threaded holes, one for location, one for impaction
- ★ Available in four sizes

No. 425 SMO Stainless Steel (B.S.S. EN 58J)
WAINWRIGHT Plate

No. 425 AL. Stainless Steel	4 in.—4 hole plate, 2½ in. blade
No. 425 A	„ 4 in.—4 hole plate, 2 in. blade
No. 425 B	„ 3 in.—3 hole plate, 2½ in. blade
No. 425 C	„ 3 in.—3 hole plate, 2 in. blade

Bilateral osteo-arthritis of the hip joint treated by McMurray Displacement Osteotomy using blade plates for internal fixation



BRIDGEND · GLAM · G.B.

Phone: BRIDGEND 938/1208 Grams: ZIMORPEDIC, BRIDGEND

London Showrooms: **Zimmer House, 176-8 Brompton Road, S.W.3**
 KNightsbridge 1919/9672