

Mol-Iron advertisement.


[s.l.]: [s.n.], 1960

<https://digital.library.wisc.edu/1711.dl/6BHT2YOJIDU6R9B>

<http://rightsstatements.org/vocab/InC/1.0/>

The libraries provide public access to a wide range of material, including online exhibits, digitized collections, archival finding aids, our catalog, online articles, and a growing range of materials in many media.

When possible, we provide rights information in catalog records, finding aids, and other metadata that accompanies collections or items. However, it is always the user's obligation to evaluate copyright and rights issues in light of their own use.



The standard by which the effectiveness
of other iron therapy **MUST** be measured

MOL-IRON[®]

a *specially processed, co-precipitated*
complex of *molybdenized* iron offering
all these important advantages:

● **MORE** hemoglobin with ● **LESS** medication
in a ● **SHORTER** period of time ● **GREATER**
patient tolerance. ● and . . . costs no more than
ordinary iron preparations.

There is a MOL-IRON product for all of your
patient needs, as listed on pp. 878 to 880 in your
1960 Physicians Desk Reference.



1 Erythrocytes 2 Polymorphonuclear Neutrophile
3 Lymphocyte 4 Monocyte 5 Eosinophile 6 Basophile
White Laboratories, Inc., Kenilworth, New Jersey