



LIBRARIES
UNIVERSITY OF WISCONSIN - MADISON

Ostic Plaster Bandages advertisement.

[s.l.]: [s.n.], 1957

<https://digital.library.wisc.edu/1711.dl/PPBQOIBB7XHVO8H>

<http://rightsstatements.org/vocab/InC/1.0/>

The libraries provide public access to a wide range of material, including online exhibits, digitized collections, archival finding aids, our catalog, online articles, and a growing range of materials in many media.

When possible, we provide rights information in catalog records, finding aids, and other metadata that accompanies collections or items. However, it is always the user's obligation to evaluate copyright and rights issues in light of their own use.

OSTIC[®] gives you something better
in a plaster bandage

THE PROOF IS IN THE ALL-DAY, EVERYDAY USE

You can feel the difference immediately
—OSTIC is creamy, gentle, never gritty

Strong, fast and precise, this is the all-around, easy-to-use bandage so popular in teaching centers. It delivers at least 97% of the original plaster to the cast because of the inseparable bonding of plaster and starch-free crinoline. Available in two settings of Fast (5 to 9 min.) and Extra-Fast (2 to 5 min.).

Ostic plaster casts stand up under all kinds of wear. Even with the more unusual fractures, as well as rambunctious little patients. And Ostic is consistent. Every time you open that waterproof package you know you're going to get the same good results you had with the last one.

For strength that meets every need, and creaminess that is kind to your hands all day, rely on OSTIC Plaster Bandages.

Curity[®]
OSTIC[®]
PLASTER BANDAGES



Ostic comes in individual airtight packages which prevent pre-setting. Plaster stays fresh indefinitely.

Bauer & Black
DIVISION OF THE KENDALL COMPANY

