

Mellaril advertisement.

[s.l.]: [s.n.], 1965

<https://digital.library.wisc.edu/1711.dl/7TY3V5MJL2SLX8H>

<http://rightsstatements.org/vocab/InC/1.0/>

The libraries provide public access to a wide range of material, including online exhibits, digitized collections, archival finding aids, our catalog, online articles, and a growing range of materials in many media.

When possible, we provide rights information in catalog records, finding aids, and other metadata that accompanies collections or items. However, it is always the user's obligation to evaluate copyright and rights issues in light of their own use.

New from Sandoz...
specifically designed for the
convenience of psychiatric-ward nurses
The Mellaril[®] Mellameter[®]
(thioridazine)

A more rapid, accurate method
of dispensing Mellaril Concentrate



Fast, Accurate—Dispenses medication to more patients in less time. One push of the plunger automatically delivers 25 mg. of Mellaril Concentrate with the same degree of accuracy achieved with careful dropper usage or assured in tablet form.

Sanitary, Safe—Completely closed system insures freedom from contamination of medication, virtually eliminates contact with work area, hands or clothing.

Easy, Convenient—Bottle can be held in one hand with index finger actuating the Mellameter. Once primed, it remains ready to use because the pump remains in bottle.

Exclusive, Disposable—Mellameter is a calibrated instrument designed to measure only Mellaril Concentrate. No need to clean pump when medication is used up as each 4 oz. bottle of concentrate comes with a new Mellameter.

For a demonstration of the new Mellameter, contact your Sandoz representative.